



Micro-Manager 1.4.22

For Compatibility: See Software Provider

Illumination System	Recognised in Micro-Manager as:	Required Micro-Manager version	USB Control	TTL Control (via external hardware)	Analog Control (via external hardware)	Further software information:	CoolLED Driver
pE-2	precisExcite	1.4 or later	Individual Wavelength control	Individual Wavelength control	N/A	<p>Click for link to Micromanager website</p> <p>pE-300 & pE-4000 device adapters are for Micro-Manager 1.4.22 only.</p> <p>For all other versions Contact CoolLED</p>	Windows XP onwards. 32-bit & 64-bit Click for link to CoolLED driver
pE-300^{white}	pE-300	1.4 or later	Individual Wavelength control	Global control	N/A		Windows 64-bit only Click for link to pE-300.dll
pE-300^{ultra}	pE-300	1.4 or later	Individual Wavelength control	Global, Individual and sequential Wavelength control	N/A		
pE-4000	pE-4000	1.4 or later	Individual Wavelength control	Individual Wavelength control	Individual Wavelength control		Windows 64-bit only Click for link to pE-4000.dll