




Screenshots for Micro-Manager 1.4.22 Imaging Software



The screenshot shows the 'Device Property Browser' window. It has a 'Refresh!' button and several checked checkboxes: 'Show cameras', 'Show shutters', 'Show stages', 'Show discrete changers', 'Show other devices', and 'Show read-only properties'. Below these is a table of properties.

Property	Value
pE4000-ChannelA	405
pE4000-ChannelB	490
pE4000-ChannelC	525
pE4000-ChannelD	635
pE4000-Description	CoolLED pE4000, Mainboard: v2.0.11 Pod: v2.0.0 Backend: v2.0.1
pE4000-Global State	0
pE4000-IntensityA	50
pE4000-IntensityB	50
pE4000-IntensityC	50
pE4000-IntensityD	50
pE4000-Lock Pod	0
pE4000-Name	pE4000
pE4000-SelectionA	0
pE4000-SelectionB	0
pE4000-SelectionC	0
pE4000-SelectionD	1

Micro-Manager pE-4000.dll



Device Property Browser

Show cameras Show read-only properties
 Show shutters
 Show stages
 Show discrete changers
 Show other devices

Refresh!

Property	Value
pE300-Description	CoolLED pE300. Mainboard: v2.0.6 Pod: v2.0.1
pE300-Global State	0
pE300-IntensityA	1
pE300-IntensityB	28
pE300-IntensityC	0
pE300-Lock Pod	0
pE300-Name	pE300
pE300-SelectionA	0
pE300-SelectionB	1
pE300-SelectionC	0

Micro-Manager pE-300.dll

PrecisExcite-IntensityA	<input checked="" type="checkbox"/>	0	<input type="text"/>
PrecisExcite-IntensityB	<input checked="" type="checkbox"/>	42	<input type="text"/>
PrecisExcite-IntensityC	<input checked="" type="checkbox"/>	0	<input type="text"/>
PrecisExcite-IntensityD	<input type="checkbox"/>	0	<input type="text"/>

precisExcite pE-2 .dll