



## Andor iQ

For Compatibility: See Software Provider

illumination System	Recognised in Andor iQ as:	Required Andor iQ version	USB Control	TTL Control (via external hardware)	Analog Control (via external hardware)	Further software information:	CoolLED Driver
<b>pE-2</b>	precisExcite	1.9.1 or later	Individual Wavelength control	Individual Wavelength control	N/A- Andor iQ offers one analog line used for piezo	<a href="#">Click for link to Andor iQ website</a>	Windows XP onwards. 32-bit & 64-bit <a href="#">Click for link to CoolLED driver</a>
<b>pE-300<sup>white</sup></b>	pE-300	3.3.2 or later	Individual Wavelength control	Global control			
<b>pE-300<sup>ultra</sup></b>	pE-300	3.3.2 or later	Individual Wavelength control	Global, Individual and sequential Wavelength control			
<b>pE-340<sup>fura</sup></b>	pE-300	3.3.2 or later	Individual Wavelength control	Global, Individual and sequential Wavelength control			
<b>pE-4000</b>	pE-4000	3.3.2 or later	Individual Wavelength control	Individual Wavelength control			