



HCImage

For Compatibility: See Software Provider

illumination System	Recognised in HCImage as:	Required HCImage version	USB Control	TTL Control (via external hardware)	Analog Control (via external hardware)	Further software information:	CoolLED Driver
pE-2	precisExcite	4.3 or later	Individual Wavelength control	Individual Wavelength control	N/A	Click for link to HCImage website	Windows XP onwards. 32-bit & 64-bit Click for link to CoolLED driver
pE-300^{white}	pE-300	4.3 or later	Individual Wavelength control	Global control	N/A		
pE-300^{ultra}	pE-300	4.3 or later	Individual Wavelength control	Global, Individual and sequential Wavelength control	N/A		
pE-340^{fura}	pE-300	4.3 or later	Individual Wavelength control	Global, Individual and sequential Wavelength control			
pE-4000	pE-4000	4.3 or later	Individual Wavelength control	Individual Wavelength control	Individual Wavelength control		