



NIS Elements

For Compatibility: See Software Provider

illumination System	Recognised in NIS Elements as:	Required NIS Elements version	USB Control	TTL Control (via external hardware)	Analog Control (via external hardware)	Further software information:	CoolLED Driver
pE-2	pE-2	4.3.02 or later	Individual Wavelength control	Individual Wavelength control	N/A	Click for link to NIS Elements website	Windows XP onwards. 32-bit & 64-bit Click for link to CoolLED driver
pE-300^{white}	pE-300	4.3.02 or later	Individual Wavelength control	Global control	N/A		
pE-300^{ultra}	pE-300	4.3.02 or later	Individual Wavelength control	Global, Individual and sequential Wavelength control	N/A		
pE-340^{fura}	pE-300	4.3.02 or later	Individual Wavelength control	Global, Individual and sequential Wavelength control	N/A		
pE-4000	pE-4000	4.3.02 or later	Individual Wavelength control	Individual Wavelength control	Individual Wavelength control		