



## Slidebook

For Compatibility: See Software Provider

Illumination System	Recognised in Slidebook as:	Required Slidebook version	USB Control	TTL Control (via external hardware)	Analogue Control (via external hardware)	Further software information:	CoolLED Driver
pE-2	precisExcite	5.0 or above	Individual Wavelength control	Individual Wavelength control	N/A	<a href="#">Click for link to Intelligent Imaging Innovations Website</a>	Windows XP onwards. 32-bit & 64-bit <a href="#">Click for link to CoolLED driver</a>
pE-300 <sup>white</sup>	pE-300	6.0.13 or above	Individual Wavelength control	Global control	N/A		
pE-300 <sup>ultra</sup>	pE-300	6.0.13 or above	Individual Wavelength control	Global, Individual and sequential Wavelength control	N/A		
pE-340 <sup>fura</sup>	pE-300	6.0.13 or above	Individual Wavelength control	Global, Individual and sequential Wavelength control	N/A		
pE-4000	pE-4000	6.0.4 or above	Individual Wavelength control	Individual Wavelength control	Individual Wavelength control		