



<b>ZEN</b>							
For Compatibility: See Software Provider							
illumination System	Recognised in ZEN as:	Required ZEN & MTB version	USB Control	TTL Control (via external hardware)	Analogue Control (via external hardware)	Further software information:	CoolLED Driver
<b>pE-300<sup>white</sup></b>	Zeiss Colibri/Colibri.2	ZEN Blue MTB 2.3.0.13 or above	Individual Wavelength control	Global control	N/A	<a href="#">Click for link to Zeiss Website</a>	Windows XP onwards. 32-bit & 64-bit <a href="#">Click for link to CoolLED driver</a>
<b>pE-300<sup>ultra</sup></b>	Zeiss Colibri/Colibri.2	ZEN Blue MTB 2.3.0.13 or above	Individual Wavelength control	Global, Individual and sequential Wavelength control	N/A		
<b>pE-340<sup>fura</sup></b>	Zeiss Colibri/Colibri.2	ZEN Blue MTB 2.3.0.13 or above	Individual Wavelength control	Global, Individual and sequential Wavelength control	N/A		
<b>pE-4000</b>	Zeiss Colibri/Colibri.2	ZEN Blue MTB 2.3.0.13 or above	Individual Wavelength control	Individual Wavelength control	Individual Wavelength control		

