



ImSpector Pro

For Compatibility: See Software Provider

Illumination System	Recognised in Slidebook as:	Required Slidebook version	USB Control	TTL Control (via external hardware)	Analogue Control (via external hardware)	Further software information:	CoolLED Driver
pE-300 ^{white}	CoolLED	5.1 or above	Individual Wavelength control	Global control	N/A	Click for link to LaVision BiotTec Website	Windows XP onwards. 32-bit & 64-bit Click for link to CoolLED driver
pE-300 ^{ultra}	CoolLED	5.1 or above	Individual Wavelength control	Global, Individual and sequential Wavelength control	N/A		
pE-340 ^{fura}	CoolLED	5.1 or above	Individual Wavelength control	Global, Individual and sequential Wavelength control	N/A		