



<h2 style="text-align: center;">NIS Elements</h2> <p style="text-align: center;">For Compatibility: See Software Provider</p>							
<b>Illumination System</b>	<b>Recognised in NIS Elements as:</b>	<b>Required NIS Elements version</b>	<b>USB Control</b>	<b>TTL Control (via external hardware)</b>	<b>Analog Control (via external hardware)</b>	<b>Further software information:</b>	<b>CoolLED Driver</b>
<b>pE-2</b>	pE-2	4.3.02 or later	Individual Wavelength control	Individual Wavelength control	N/A	<a href="#">Click for link to NIS Elements website</a>	Windows XP onwards. 32-bit & 64-bit <a href="#">Click for link to CoolLED driver</a>
<b>pE-300<sup>white</sup></b>	pE-300	4.40 or later	Individual Wavelength control	Global control	N/A		
<b>pE-300<sup>ultra</sup></b>	pE-300	4.40 or later	Individual Wavelength control	Global, Individual and sequential Wavelength control	N/A		
<b>pE-340<sup>fura</sup></b>	pE-300	4.40 or later	Individual Wavelength control	Global, Individual and sequential Wavelength control	N/A		
<b>pE-4000</b>	pE-4000	4.40 or later	Individual Wavelength control	Individual Wavelength control	Individual Wavelength control		
<b>pE-800</b>	pE-800	v5.30.01 + pE-800 patch or v5.30.02	Individual Wavelength control	Global, Individual and sequential Wavelength control	Individual Wavelength control	<a href="#">NIS pE-800 Patch</a>  <a href="#">pE-800 Patch Install Guides</a>	Windows 7 onwards. Pre-Windows 10 requires <a href="#">CoolLED driver</a>