



pE-340^{fura}

pE-340^{fura}
pE-340-FR-D-YYY-ZZ



pE-340^{fura} light guide
pE-340-FR-L-ZZ



Key :

- 1 Power supply
- 2 Control Pod
- 3 pE-340^{fura} Direct Fit
- 4 pE-340^{fura} light guide
- 5 pE-Adaptor
- 6 3 m liquid light guide
- 7 pE- Universal Collimator
- 8 Global TTL
- 9 Individual channel TTL
- 10 USB
11. Removable inline excitation filter holders