



HCIImage							
For Compatibility: See Software Provider							
illumination System/Accessory	Recognised in HCIImage as:	Required HCIImage version	USB Control	TTL Control (via external hardware)	Analog Control (via external hardware)	Further software information:	CoolLED Driver
pE-2	precisExcite	4.0 or later	Individual Wavelength control	Individual Wavelength control	N/A	Click for link to HCIImage website	Windows XP onwards. 32-bit & 64-bit Click for link to CoolLED driver
pE-300^{white}	pE-300	4.0 or later	Individual Wavelength control	Global control	N/A		
pE-300^{ultra}	pE-300	4.0 or later	Individual Wavelength control	Global, Individual and sequential Wavelength control	N/A		
pE-340^{fura}	pE-300	4.0 or later	Individual Wavelength control	Global, Individual and sequential Wavelength control	N/A		
pE-4000	pE-4000	4.4 or later	Individual Wavelength control	Individual and Wavelength control	Individual Wavelength control		
pE-800	pE-800	4.7 or later	Individual Wavelength control	Global, Individual and sequential Wavelength control	Individual Wavelength control		
pE-6501-8	NI 6501	4.7 or later	N/A	N/A	N/A		