



## HCImage

For Compatibility: See Software Provider

illumination System/Accessory	Recognised in HCImage as:	Required HCImage version	USB Control	TTL Control (via external hardware)	Analog Control (via external hardware)	Further software information:	CoolLED Driver
pE-2	precisExcite	4.0 or later	Individual Wavelength control	Individual Wavelength control	N/A	<a href="#">Click for link to HCImage website</a>	Windows XP onwards. 32-bit & 64-bit <a href="#">Click for link to CoolLED driver</a>
pE-300 <sup>white</sup>	pE-300	4.0 or later	Individual Wavelength control	Global control	N/A		
pE-300 <sup>ultra</sup>	pE-300	4.0 or later	Individual Wavelength control	Global, Individual and sequential Wavelength control	N/A		
pE-340 <sup>fura</sup>	pE-300	4.0 or later	Individual Wavelength control	Global, Individual and sequential Wavelength control	N/A		
pE-4000	pE-4000	4.4 or later	Individual Wavelength control	Individual and Wavelength control	Individual Wavelength control		
pE-800 Series	pE-800	4.7 or later	Individual Wavelength control	Global, Individual and sequential Wavelength control	Individual Wavelength control		
pE-6501 pE-6501-8 pE-6501-8-FR	NI 6501	4.7 or later	N/A	N/A	N/A		