



## ImSpector Pro

For Compatibility: See Software Provider

illumination System	Recognised in ImSpector Pro as:	Required ImSpector Pro version	USB Control	TTL Control (via external hardware)	Analogue Control (via external hardware)	Further software information:	CoolLED Driver
pE-300 <sup>white</sup>	CoolLED	5.1 or above	Individual Wavelength control	Global control	N/A	<a href="#">Click for link to Miltenyi BioTec Website</a>	Windows XP onwards. 32-bit & 64-bit <a href="#">Click for link to CoolLED driver</a>
pE-300 <sup>ultra</sup>	CoolLED	5.1 or above	Individual Wavelength control	Global, Individual and sequential Wavelength control	N/A		
pE-340 <sup>fura</sup>	CoolLED	5.1 or above	Individual Wavelength control	Global, Individual and sequential Wavelength control	N/A		