



## Inscoper

For Compatibility: See Software Provider

illumination System	Recognised in INSCOPER as:	Required INSCOPER version	USB Control	TTL Control (natively integrated)	Analog Control (natively integrated)	Further software information:	CoolLED Driver
pE-2	pE-2	From version 1.0 onwards	Individual Wavelength control	Individual Wavelength control	N/A	<a href="#">Click for link to Inscoper website</a>	Not required
pE-300 <sup>white</sup>	pE-300white		Individual Wavelength control	Global control	N/A		
pE-300 <sup>ultra</sup>	pE-300ultra		Individual Wavelength control	Individual Wavelength control	N/A		
pE-340 <sup>fura</sup>	pE-340fura		Individual Wavelength control	Individual Wavelength control	N/A		
pE-800 Series	pE-800		Individual Wavelength control	Individual Wavelength control	Individual Wavelength control		
pE-4000	pE-4000		Individual Wavelength control	Individual Wavelength control	Individual Wavelength control		