



MetaMorph

For Compatibility: See Software Provider

illumination System	Recognised in MetaMorph as:	Required MetaMorph version	USB Control	TTL Control (via external hardware)	Analog Control (via external hardware)	Further software information:	CoolLED Driver
pE-2	precisExcite	Windows XP: 7.1.3 Windows 7 onwards: 7.7.2	Individual Wavelength control	Individual Wavelength control	N/A	Click for link to MolDev website N.B. We have experienced compatibility issues with the following Meta Imaging versions: 7.8.10 7.8.11 7.8.12	Windows XP onwards. 32-bit & 64-bit Click for link to CoolLED driver
pE-300 ^{white}	pE-4000	7.8.9 or later	Individual Wavelength control	Global control	N/A		
pE-300 ^{ultra}	pE-4000	7.8.9 or later	Individual Wavelength control	Global, Individual and sequential Wavelength control	N/A		
pE-340 ^{fura}	pE-4000	7.8.9 or later	Individual Wavelength control	Global, Individual and sequential Wavelength control	N/A		
pE-4000	pE-4000	7.8.9 or later	Individual Wavelength control	Individual Wavelength control	Individual Wavelength control		
pE-6501-8	National Instruments NI-DAQmx	7.8.9 or later	N/A	N/A	N/A		