



<h2 style="text-align: center;">MetaMorph</h2> <p style="text-align: center;">For Compatibility: See Software Provider</p>							
Illumination System	Recognised in MetaMorph as:	Required MetaMorph version	USB Control	TTL Control (via external hardware)	Analog Control (via external hardware)	Further software information:	CoolLED Driver
pE-2	CoolLED pE	Windows XP: 7.1.3 Windows 7 onwards: 7.7.2	Individual Wavelength control	Individual Wavelength control	N/A	<p>N.B. We have experienced compatibility issues with the following Meta Imaging versions:</p> <p>7.8.10 7.8.11 7.8.12</p>	<p>Windows XP onwards. 32-bit &amp; 64-bit</p> <p><a href="#">Click for link to CoolLED driver</a></p>
pE-300 <sup>white</sup>	CoolLED pE	7.8.9 or later	Individual Wavelength control	Global control	N/A		
pE-300 <sup>ultra</sup>	CoolLED pE	7.8.9 or later	Individual Wavelength control	Global, Individual and sequential Wavelength control	N/A		
pE-340 <sup>fura</sup>	CoolLED pE	7.8.9 or later	Individual Wavelength control	Global, Individual and sequential Wavelength control	N/A		
pE-4000	CoolLED pE	7.8.9 or later	Individual Wavelength control	Individual Wavelength control	Individual Wavelength control		
pE-800 Series	CoolLED pE	32-bit not supported Windows 10 onwards: 7.10.4 ++	Individual Wavelength control	Global, Individual and sequential Wavelength control	Individual Wavelength control		
pE-6501 pE-6501-8 pE-6501-8-FR	National Instruments NI-DAQmx	7.8.9 or later	N/A	N/A	N/A		