



| <h2 style="text-align: center;">NIS Elements</h2> <p style="text-align: center;">For Compatibility: See Software Provider</p> | | | | | | | |
|---|---------------------------------------|--------------------------------------|-------------------------------|--|---|---|---|
| Illumination System | Recognised in NIS Elements as: | Required NIS Elements version | USB Control | TTL Control (via external hardware) | Analog Control (via external hardware) | Further software information: | CoolLED Driver |
| pE-2 | pE-2 | 4.3.02 or later | Individual Wavelength control | Individual Wavelength control | N/A | Click for link to NIS Elements website | Windows XP onwards. 32-bit & 64-bit Click for link to CoolLED driver |
| pE-300^{white} | pE-300 | 4.40 or later | Individual Wavelength control | Global control | N/A | | |
| pE-300^{ultra} | pE-300 | 4.40 or later | Individual Wavelength control | Global, Individual and sequential Wavelength control | N/A | | |
| pE-340^{fura} | pE-300 | 4.40 or later | Individual Wavelength control | Global, Individual and sequential Wavelength control | N/A | | |
| pE-4000 | pE-4000 | 4.40 or later | Individual Wavelength control | Individual Wavelength control | Individual Wavelength control | | |
| pE-800 Series (pE-800 & pE-800^{fura}) | pE-800 | v5.30.01 + pE-800 patch or v5.30.02 | Individual Wavelength control | Global, Individual and sequential Wavelength control | Individual Wavelength control | NIS pE-800 Patch pE-800 Patch Install Guides | Windows 7 onwards. Pre-Windows 10 requires CoolLED driver |