



Slidebook

For Compatibility: See Software Provider

illumination System	Recognised in Slidebook as:	Required Slidebook version	USB Control	TTL Control (via external hardware)	Analogue Control (via external hardware)	Further software information:	CoolLED Driver
pE-2	precisExcite	5.0 or above	Individual Wavelength control	Individual Wavelength control	N/A	Click for link to Intelligent Imaging Innovations Website	Windows XP onwards. 32-bit & 64-bit Click for link to CoolLED driver
pE-300^{white}	pE-300	6.0.13 or above	Individual Wavelength control	Global control	N/A		
pE-300^{ultra}	pE-300	6.0.13 or above	Individual Wavelength control	Global, Individual and sequential Wavelength control	N/A		
pE-340^{fura}	pE-300	6.0.13 or above	Individual Wavelength control	Global, Individual and sequential Wavelength control	N/A		
pE-4000	pE-4000	6.0.4 or above	Individual Wavelength control	Individual Wavelength control	Individual Wavelength control		
pE-800 Series	pE-800	6.0.23 onwards	Individual Wavelength control	Global, Individual and sequential Wavelength control	Individual Wavelength control		