

## μManager

For Compatibility: See Software Provider

Illumination System	Recognised in μManager as:	Required μManager version	USB Control	TTL Control (via external hardware)	Analog Control (via external hardware)	Further software information:	CoolLED Driver
pE-2	precisExcite	1.4 or later	Individual Wavelength control	Individual Wavelength control	N/A	<a href="#">Click for link to Micromanager website</a>  pE-300 & pE-4000 device adapters are present in μManager <b>1.4.23</b> nightly builds as standard  The device adapters can also be downloaded and added to the official <b>1.4.22</b> release only.	Windows XP onwards. 32-bit & 64-bit <a href="#">Click for link to CoolLED driver</a>
pE-300 <sup>white</sup>	pE-300	1.4 or later	Individual Wavelength control	Global control	N/A		Windows 64-bit only <a href="#">Click for link to pE-300.dll</a>
pE-300 <sup>ultra</sup>	pE-300	1.4 or later	Individual Wavelength control	Global, Individual and sequential Wavelength control	N/A		Windows 64-bit only <a href="#">Click for link to pE-4000.dll</a>
pE-4000	pE-4000	1.4 or later	Individual Wavelength control	Individual Wavelength control	Individual Wavelength control		Windows 64-bit only <a href="#">Click for link to pE-800.dll</a>
pE-800 Series	pE-800	1.4.22	Individual Wavelength control	Individual Wavelength control	Individual Wavelength control		Windows 64-bit only <a href="#">Click for link to pE-300.dll</a>
pE-340 <sup>fura</sup>	pE-300	1.4 or later	Individual Wavelength control	Global, Individual and sequential Wavelength control	N/A		N/A
pE-6501 pE-6501-8 pE-6501-8-FR	NI100X	1.4 or later	N/A	N/A	N/A		<a href="#">pE-6501 Guide</a> <a href="#">pE-6501-8 Guide</a>