

μManager

For Compatibility: See Software Provider

illumination System	Recognised in μManager as:	Required μManager version	USB Control	TTL Control (via external hardware)	Analog Control (via external hardware)	Further software information:	CoolLED Driver
pE-2	precisExcite	1.4 or later	Individual Wavelength control	Individual Wavelength control	N/A	Click for link to Micromanager website	Windows XP onwards. 32-bit & 64-bit Click for link to CoolLED driver
pE-300 ^{white}	pE-300	1.4 or later	Individual Wavelength control	Global control	N/A		pE-300 & pE-4000 device adapters are present in μManager 1.4.23 nightly builds as standard
pE-300 ^{ultra}	pE-300	1.4 or later	Individual Wavelength control	Global, Individual and sequential Wavelength control	N/A	The device adapters can also be downloaded and added to the official 1.4.22 release only.	Windows 64-bit only Click for link to pE-4000.dll
pE-4000	pE-4000	1.4 or later	Individual Wavelength control	Individual Wavelength control	Individual Wavelength control		Windows 64-bit only Click for link to pE-800.dll
pE-800 Series	pE-800	1.4.22	Individual Wavelength control	Individual Wavelength control	Individual Wavelength control		N/A
pE-340 ^{fura}	Not recognised	N/A	N/A	N/A	N/A	N/A	N/A
pE-6501 pE-6501-8 pE-6501-8-FR	NI100X	1.4 or later	N/A	N/A	N/A	pE-6501 Guide pE-6501-8 Guide	N/A