



Volocity

For Compatibility: See Software Provider

Illumination System/Accessory	Recognised in Volocity as:	Required Volocity version	USB Control	TTL Control (via external hardware)	Analog Control (via external hardware)	Further software information:	CoolLED Driver
pE-300 ^{white}	CoolLED	6.5.1 or later	Individual Wavelength control	Global control	N/A	Click here for Quorum Technologies Volocity website	Windows XP onwards. 32-bit & 64-bit Click for link to CoolLED driver
pE-300 ^{ultra}	CoolLED	6.5.1 or later	Individual Wavelength control	Global, Individual and sequential Wavelength control	N/A		
pE-340 ^{fura}	CoolLED	6.5.1 or later	Individual Wavelength control	Global, Individual and sequential Wavelength control	N/A		
pE-4000	CoolLED	6.5.1 or later	Individual Wavelength control	Individual Wavelength control	Individual Wavelength control		