

# pT-100

pT-100  
pT-100-XXX-YYY-Z



pT-100 with 1.5 m Liquid Light Guide  
pT-100-XXX-L15-ZZ



pT-100 with 3 m Liquid Light Guide  
pT-100-XXX-L3-ZZ



Transmitted Port

**Key :**

- 1: Power supply
- 2: Control Pod
- 3: pT-100
- 4: pT-100 Liquid Light Guide 1.5 m long
- 5: pE-Universal Collimator
- 6: pE-Adaptor
- 7: pT-100 Liquid Light Guide 3 m long