



# pE-400

pE-400 direct fit  
pE-400-D-SB-YYY-ZZ  
pE-400-D-MB-YYY-ZZ  
(SB = 365 nm, MB = 365 nm, MB = 400 nm)



pE-400 light guide  
pE-400-L-SB-SYS-ZZ  
pE-400-L-MB-SYS-ZZ  
(SB = 365 nm, MB = 400 nm)



### Key :

- 1 Power supply
- 2 Control Pod
- 3 pE-400 direct fit
- 4 pE-400 light guide
- 5 pE-Adaptor
- 6 1.5 m liquid light guide
- 7 3 m liquid light guide
- 8 pE-Universal Collimator
- 9 Global TTL
- 10 USB