



## Media Cybernetics Image-Pro

Illumination System	Recognised in Image-Pro as	Required Image-Pro version	USB Control	TTL Control (via external hardware)	Analogue Control (via external hardware)	CoolLED Driver	Further software information
<b>pE-300<sup>ultra</sup></b>	pE-300	Image-Pro v11	Individual channel control	No	N/A	Windows XP onwards. 32-bit & 64-bit All Windows versions require <a href="#">CoolLED driver</a>	<a href="#">Image-Pro website</a>  <a href="#">View Quick Start Guide</a>
<b>pE-300<sup>white</sup></b>	pE-300	Image-Pro v11	Individual channel control	No	N/A		
<b>pE-340<sup>fura</sup></b>	pE-300	Image-Pro v11	Individual channel control	No	N/A		