



Media Cybernetics Image-Pro For Compatibility: See Software Provider

Illumination System	Recognised in Image-Pro as:	Required Image-Pro version	USB Control	TTL Control (via external hardware)	Analogue Control (via external hardware)	CoolLED Driver	Further software information:
pE-300^{ultra}	pE-300	Image-Pro v11	Individual channel control	No	N/A	Windows XP onwards. 32-bit & 64-bit Driver download	Visit Image-Pro website
pE-300^{white}	pE-300	Image-Pro v11	Individual channel control	No	N/A		View Quick Start Guide
pE-340^{fura}	pE-300	Image-Pro v11	Individual channel control	No	N/A		